


ENFit[®] Enteral Syringes

Comply with ISO80369-3 standard



The worldwide introduction of the new ENFit connector will reduce the risk of enteral tube feeding misconnections and improve patient safety throughout hospitals around the globe.

Usage of ENFit Enteral Syringes

• Enteral syringes can be used for delivery of enteral feeding, aspiration, flushing and to deliver enteral medicines to all patients groups requiring assistance for feeding, whether using a standard feeding bag, formula bottle or syringe pump.



◆ Sizes available:

1ml low dose tip, 2.5ml low dose tip, 5ml low dose tip, 10ml, 20ml, 30ml, 60ml and 100ml

Advantage of our ENFit Enteral Syringes

- Super quality by using medical grade PP
- Transparency of barrel
- Strong ink adhesion to ensure legible and clear graduation
- Latex-free piston
- Using silicone oil of medical grade
- Free of pyrogen and hemolysis
- DEHP free
- Strict quality control by QC department in every process
- Fully compliant with ISO standard
- OEM accepted



*** Reusable ENFit Enteral Syringes and Oral Tip Syringes are also available.**

Shantou Wealy Medical Instrument Co., Ltd.

North Jinhuan Road (near of Qishan Mid-school) 515064 Shantou, China

Tel: +86-754-82121967

Fax: +86-754-82121654

Website: www.wealy.com

E-mail: export@wealy.com

Wealy[®]